beurer

BF 600 Style Pure black & white

DE Diagnosewaage Gebrauchsanweisung	2	IT Bilancia diagnostica Istruzioni per l'uso	53
EN Diagnostic bathroom scale Instructions for use	15	TR Diyagnoz terazisi Kullanım kılavuzu	65
FR Pèse-personne impédancemètre Mode d'emploi	27	RU Диагностические весы Инструкция по применению	77
ES Báscula de diagnóstico Instrucciones de uso	40	PL Waga diagnostyczna Instrukcja obsługi	91

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www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php

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Read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Dear customer,

Thank you for choosing a product from our range. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty, baby and air.

With kind regards,
Your Beurer team

Contents

1. Getting to know your device	15	9. Take measurement	20
2. Signs and symbols	16	10. Evaluating the results	21
3. Intended use.....	17	11. Other functions.....	23
4. Warnings and safety notes	17	12. Cleaning and maintaining the device	24
5. Information	18	13. What if there are problems?	24
6. Device description.....	19	14. Cleaning and maintenance.....	25
7. Initial use with the app.....	19	15. Technical specifications	25
8. Initial use of the scale without the app.....	20	16. Warranty/service.....	26

Included in delivery

- Quick guide
- BF 600 diagnostic bathroom scale
- 2 x 1.5 V AAA batteries
- These instructions for use
- 1x warranty leaflet

1. Getting to know your device

Functions of the device

This digital diagnostic scale is intended for weighing and providing a diagnosis of your personal fitness data.

It is intended for self-testing in the private domain.

The scale is equipped with the following diagnostic functions, which can be used by up to 8 different people:

- Body weight measurement and BMI,
- Calculating body fat percentage,
- Body water content,
- Muscle percentage,
- Bone mass and
- Basal metabolic rate and active metabolic rate.

The scale also includes the following additional features:

- Switch between kilograms “kg”, pounds “lb” and stones “st”,
- Automatic switch-off function,
- Low battery indicator,
- Automatic user recognition,
- Storage of the last 30 measurements for 8 users if the data cannot be transferred to the app.

System requirements „beurer HealthManager Pro“ / „beurer HealthManager“ App



iOS from version 12.0, Android™ from version 8.0 with Bluetooth® ≥ 4.0.

List of compatible devices:

2. Signs and symbols

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	WARNING Warning notice indicating a risk of injury or damage to health.		
	IMPORTANT Safety note indicating possible damage to the device/accessory.		
	Note Note on important information		
	Observe the instructions for use		
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE		
	Dispose of packaging in an environmentally friendly manner		
	This product satisfies the requirements of the applicable European and national directives.		
	United Kingdom Conformity Assessed Mark		
	The products demonstrably meet the requirements of the Technical Regulations of the EAEU.		
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.		
	Manufacturer		
	Do not load the scale beyond 180 kg/396 lb/28 st.	 	Danger of slipping: do not step on the scale with wet feet.
 	Danger of tipping: position yourself in the centre of the weighing surface.		Place the scale on an even surface. No carpet.

3. Intended use

The device is only intended for weighing humans and for recording your personal fitness data. The device is only intended for private use, and not for medical or commercial purposes.

4. Warnings and safety notes

Please read these instructions for use carefully and keep them for later use, be sure to also make them accessible to other users and observe the notes they contain.

WARNING

- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- Do not use during pregnancy.
- Never step onto the outer edge of the scale with one foot as this could cause the scale to tip!
- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a battery be swallowed, seek medical assistance immediately.
- Keep packaging material away from children (risk of suffocation).
- Warning, do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- **⚠ Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- **⚠ Risk of explosion!** Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

General notes

- The device is suitable only for private use, and not for medical or commercial purposes.
- Please note that measuring tolerances are possible for technical reasons as this scale is not calibrated for use in a professional medical context.
- In cases of differing measurements (between the scale and app), only use the measured values displayed on the scale.
- The scale's maximum capacity is 180 kg (396 lb/28 st). The results for weight and bone mass are displayed in 100 g increments (0.2 lb).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- Calorie requirement is indicated in steps of 1 kcal.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm". You can change the unit settings during initial use of the app.
- Place the scale on an even, hard surface; a hard surface is crucial to accurate measurements.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Repairs must only be carried out by Beurer Customer Services or authorised retailers. Before making a claim, please check the batteries first and replace them if necessary.

Storage and maintenance

The accuracy of the measurements and service life of the device depend on its careful handling:

IMPORTANT

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when not in use.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Do not press the buttons violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile telephones).

5. Information

The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content within seconds, using a current that is not noticeable, is completely non-hazardous and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

Please note that the values calculated by the diagnostic scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

General tips

- As far as possible, weigh yourself at the same time everyday (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be made when bare-foot; the soles of the feet may be lightly dampened where appropriate.
Results may be unsatisfactory if the soles of the feet are completely dry or feature high levels of horny skin as this results in insufficient conductivity.

Note: There must be no skin contact between your feet, legs, calves or thighs. Otherwise, the measurement cannot be taken correctly.

- Stand up straight and still whilst weighing yourself.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in the body can disperse.
- You must remember that only the long-term trend is important. Brief deviations in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.

Limitations


When calculating the body fat and the additional values, deviating and implausible results may occur in the case of:

- Children under 10 years of age,
- Performance athletes and body-builders,
- Persons with a fever, persons undergoing dialysis, persons demonstrating oedema symptoms and persons suffering from osteoporosis,
- Persons taking cardiovascular medicine (affecting the heart and vascular system),
- Persons taking vasodilative or vasoconstrictive medication,
- Persons with significant anatomical deviations in their legs compared with their overall body height (leg length significantly shorter or longer).

6. Device description

Rear

1. Battery compartment
2. Unit/reset button

Display

3. Bluetooth® symbol for connection between scale and smartphone.
4. Measurements

7. Initial use with the app

To calculate your body fat percentage and additional body values, you must save the personal user data on the scale.

The scale has 8 user memory spaces which allows, for example, you and your family members to save personal settings.

It is also possible to activate users via other mobile end devices on which the “beurer HealthManager Pro“ / “beurer HealthManager“ app has been installed or by changing the user profile in the app (see the app settings).

Keep the smartphone close to the scale in order to maintain an active *Bluetooth*® connection during initial use.

– **Activate *Bluetooth*® in your smartphone’s settings.**

– **Insert the batteries and place the scale on a solid surface (“App” will flash on the scale’s LCD).**

– **Download the free “beurer HealthManager Pro“ / “beurer HealthManager“ app from the Apple App Store (iOS) or at Google Play (Android™).**

– **Start the app and follow the instructions.**

– **Enter the SR BF 1 related data.**

The following settings must be set or entered in the “beurer HealthManager Pro“ / “beurer HealthManager“ app:

User data	Set values
Gender	Male (♂), female (♀)
Height	100 to 220 cm (3' 3.5" to 7' 2.5")
Age	10 to 100 years
Activity level	1 to 5

Activity levels

The medium-term and long-term perspective is decisive when selecting the activity level.

Activity level	Physical activity
1	None.
2	Low: little and light physical exertion (e.g. going for walks, light gardening work, gymnastic exercises).
3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each session.
4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each session.
5	Very high: intensive physical exertion, intensive training or hard physical work, daily for at least 1 hour.

The app assigns the next free scale user memory (P-01 – P-08) to you.

– Taking the measurement

For automatic personal recognition, the first measurement must be assigned to your personal user data. The app will prompt you to do this. Step onto the scale, ideally with bare feet, and ensure that you are standing still with equal weight distribution and with both feet on the electrodes. Now continue with chapter “Take measurement”.

8. Initial use of the scale without the app

We recommend using the app for the initial use and implementing all settings. See chapter 6 “Initial use with the app” However, this can also be carried out on the device itself.

Remove the batteries from the protective packaging and insert the batteries into the scale ensuring that the battery polarity is correct. If the scale does not function, remove the batteries completely and reinsert them.

Place the scale on an even, hard surface; a hard surface is crucial to accurate measurements.

Press the “SET” button.

Set the year/month/day and hour/minute by pressing the buttons “ \vee \wedge ” and confirm these settings by pressing the “SET” button.

The scale displays “0.0”. Press the “SET” button again. “P-01” or the last used scale user memory will appear. Select your user memory (P-01 to P-08) with the arrows “ \vee \wedge ” and press the “SET” button. Carry out your settings, such as gender, height etc. (see chapter “Initial use with the app (table)”). Confirm each setting by pressing the “SET” button. A PIN will be displayed at the end of your settings. This will be required for logging into the app later. Afterwards “0.0” is displayed.

For automatic personal recognition, the first measurement must be assigned to your personal user data. Step onto the scale, ideally with bare feet, and ensure that you are standing still with equal weight distribution and with both feet on the electrodes.

7.1 Logging in to the app at a later time




If you have carried out the initial use of the scale without using the app, you can transfer or create your user on the app as follows.

– Download the free “beurer HealthManager Pro” / “beurer HealthManager” app from the Apple App Store (iOS) or at Google Play (Android™).

– Start the app and follow the instructions.

9. Take measurement

Place the scale on an even, hard surface; a hard surface is crucial to accurate measurements.

9.1 Diagnostic functions

Step onto the scale **with bare feet** and ensure that you are standing still with equal weight distribution and with both legs on the electrodes.

 **Note:** The measurement result will be incorrect if the measurement is taken with socks on.

The scale begins to take measurements straight away.

First, the weight is displayed.

During the BIA analysis, “0000” is displayed.

If a user was assigned, P-01–P-08 is displayed and all of the following values will be saved.

1. BMI
2. Body fat in %
3. Body water in % \approx
4. Muscle percentage in % \rightleftharpoons
5. Bone mineral amount in kg \rightleftharpoons
6. Basal metabolic rate in kcal (BMR)
7. Active metabolic rate in kcal (AMR)

i Note: If the user is not recognised, then the user assignment (e.g. P-01) is missing and only the weight can be displayed. The measured values will not be saved. This happens if your measured weight differs from the previous measurement by ± 3 kg or if another user has similar values (± 3 kg) to yours. In this case, repeat the measurement by turning on the scale and pressing the “^” button to select your user memory (P-01–P-08). Repeat the measurement.

9.2 Weight

Step onto the scale **wearing shoes** and ensure that you are standing still with equal weight distribution and with both legs on the scale. The scale begins to take measurements straight away.

First, the weight is displayed.

“0” will briefly be displayed on the LCD. A BIA analysis is not possible.

If a user was assigned, “P-01”, for example, will be displayed and the weight will be saved.

i Note: If the user is not recognised, then the user assignment (e.g. P-01) is missing and only the weight can be displayed. The measured values will not be saved. This happens if your measured weight differs from the previous measurement by ± 3 kg or if another user has similar values (± 3 kg) to yours. In this case, repeat the measurement by turning on the scale and pressing the “^” button to select your user memory (P-01–P-08). Repeat the measurement.

10. Evaluating the results

Body mass index

The Body Mass Index (BMI) is a figure that is commonly used to evaluate body weight. The figure is calculated from body weight and height. The formula is:

Body mass index = body weight \div height². The measurement unit for BMI is [kg/m²]. According to the BMI, weight is classified for adults (20 years and over) using the following values:

Category		BMI
Underweight	Severely underweight	<16
	Underweight	16-16.9
	Slightly underweight	17-18.4
Normal weight		18.5-25
Overweight	Overweight	25.1-29.9
Obese (overweight)	Class I obesity	30-34.9
	Class II obesity	35-39.9
	Class III obesity	≥ 40

Body fat percentage

The subsequent body fat percentages serve as a guideline (for further information ask your doctor).

Male

Age	Low	Normal	Moderate	High
10-14	<11%	11-16%	16.1-21%	>21.1%
15-19	<12%	12-17%	17.1-22%	>22.1%
20-29	<13%	13-18%	18.1-23%	>23.1%
30-39	<14%	14-19%	19.1-24%	>24.1%
40-49	<15%	15-20%	20.1-25%	>25.1%
50-59	<16%	16-21%	21.1-26%	>26.1%
60-69	<17%	17-22%	22.1-27%	>27.1%
70-100	<18%	18-23%	23.1-28%	>28.1%

Female

Age	Low	Normal	Moderate	High
10-14	<16%	16-21%	21.1-26%	>26.1%
15-19	<17%	17-22%	22.1-27%	>27.1%
20-29	<18%	18-23%	23.1-28%	>28.1%
30-39	<19%	19-24%	24.1-29%	>29.1%
40-49	<20%	20-25%	25.1-30%	>30.1%
50-59	<21%	21-26%	26.1-31%	>31.1%
60-69	<22%	22-27%	27.1-32%	>32.1%
70-100	<23%	23-28%	28.1-33%	>33.1%

A lower value will often be determined in the case of athletes. Depending on the sporting activities performed, the training intensity and the physical constitution, the results may still be below the stated standard values. Please note, however, that there may be a risk to health in the case of extremely low values.

Body water content

The body water content in % is normally within the following ranges:

Male

Age	Poor	Good	Very good
10-100	<50 %	50-65%	>65%

Female

Age	Poor	Good	Very good
10-100	<45 %	45-60%	>60%

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle percentage

The muscle percentage is normally within the following ranges:

Male

Age	Low	Normal	High
10-14	<44 %	44-57%	>57%
15-19	<43 %	43-56%	>56%
20-29	<42%	42-54%	>54%
30-39	<41%	41-52%	>52%
40-49	<40%	40-50%	>50%
50-59	<39%	39-48%	>48%
60-69	<38%	38-47%	>47%
70-100	<37 %	37-46%	>46%

Female

Age	Low	Normal	High
10-14	<36 %	36-43%	>43%
15-19	<35 %	35-41%	>41%
20-29	<34%	34-39%	>39%
30-39	<33%	33-38%	>38%
40-49	<31%	31-36%	>36%
50-59	<29%	29-34%	>34%
60-69	<28%	28-33%	>33%
70-100	<27 %	27-32%	>32%

Bone mass

Like the rest of our body, our bones undergo natural processes of growth, breakdown and ageing. Bone mass increases rapidly during childhood and reaches its peak between 30 and 40 years. As you age, your bone mass starts to gradually decline. If you have a healthy diet (particularly calcium and vitamin D) and regularly exercise, you can slow this process down. You can also strengthen the stability of your skeletal structure using targeted muscle growth. Please note that this scale does measure minerals of bones. Bone mass is very difficult to modify but does vary slightly as a result of influencing factors (weight, height, age, sex). There are no recognised guidelines and recommendations available.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. when lying in bed 24 hours a day). This value is primarily dependent on weight, height and age.

It is displayed on the diagnostic scale in the unit kcal/day and calculated using the scientifically recognised Harris-Benedict equation.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic scale calculates this by means of the activity level (1–5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in than is used over a long period of time, the body takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Temporal context of the results

 Please note that only the long-term trend is important. Brief deviations in weight within a few days are normally caused by loss of fluids.

The interpretation of the results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as on the length of time over which these changes occur.

Rapid changes within the scope of a few days are to be considered separately from medium-term changes (in the scope of weeks) and long-term changes (months).

It can be said as a basic rule that short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may also relate to the fat percentage and the muscle percentage.

- If your weight drops in the short term but your body fat percentage increases or stays the same, you have only lost water, for example following a training session, visit to the sauna or a diet aimed only at fast weight loss.
- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have built up valuable muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working – you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water and muscle percentages must not be added (muscle tissue also contains components made of body water).

11. Other functions

User assignment

It is possible to assign up to 8 scale users for a weight-only measurement (with shoes) and for a diagnosis measurement (bare feet).

For a new measurement, the scale assigns the measurement to the user whose most recently saved measurement falls within $\pm 3\text{kg}$ of this. See chapter “Take measurement”.

Storing the assigned measurements on the scale and transferring them to the app.

If the app is not open, the newly assigned measurements are saved on the scale. A total of 30 measurements per user can be saved on the scale. The saved measurements are transferred automatically to the app if you open the app within the *Bluetooth*[®] range.

It is possible to automatically transfer data when the scale is switched off.

Deleting scale data

If you would like to completely delete all measurements and old user data on the scale, remove the batteries for at least 10 seconds. Then press the unit/reset button and reinsert the batteries at the same time. “.” will appear on the display for several seconds.

This is necessary if you would like to delete scale users who have been set up incorrectly or are no longer required, for example.

You must then carry out the actions in chapter “Initial use with the app” again (user assignment).

Changing the batteries

Your scale is equipped with a low battery indicator. If you operate the scale with flat batteries, “Lo” will appear on the display and the scale will automatically switch off. In this case, the batteries must be replaced (2 x 1.5 V AAA). All saved measurements and users remain saved.

Connect the scale to the app again. The system time is updated.

NOTE:

- When changing the batteries, use batteries of the same type, make and capacity.
- Do not use rechargeable batteries.
- Use batteries free from heavy metals.

12. Cleaning and maintaining the device

The device should be cleaned from time to time.

This should be done using a damp cloth and, if necessary, a small amount of detergent.

IMPORTANT

- Never use abrasive solvents or cleaning products!
- Never submerge the device in water!
- Do not clean the device in a dishwasher!

13. What if there are problems?

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
Only weight display and P-01 - P-08 are missing	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Select user on the scale using the “^” button or repeat user assignment in the app. See chapter “Other functions”.
Err	The proportion of fat is outside the measurable range (less than 4 % or greater than 65%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary. If Err is still displayed, you are outside of the measurable range.
Old	Maximum weight capacity of 180 kg exceeded.	Loads must be limited to 180 kg.
Incorrect weight is displayed.	No flat, stable surface.	Place the scale on a flat, stable surface. No carpet.
	Not standing still.	Stand as still as possible.
	Incorrect zero point of scale.	Wait until the scale switches itself off. Activate scale. “C” may be displayed. Wait for “0.0 kg”. Then repeat the measurement.

Display	Cause	Solution
No Bluetooth ® connection (📶 symbol is missing).	Device outside the range.	The range in the open is approx. 15 m. Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as WIFI routers, microwaves, induction hobs etc.
	No connection to the app.	Completely close the app (including in the background). Switch Bluetooth ® off and back on again. Switch the smartphone off and back on again. Remove the batteries in the scale for a short period and reinsert them. See FAQs at www.beurer.com .
FULL	User memory space is full. New measurements overwrite the older measurements.	Open the app. The data is transferred automatically. This may take up to one minute.
Lo	The batteries in the scale are empty.	Replace the batteries in the scale.

14. Cleaning and maintenance

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

For environmental reasons, do not dispose of the device in the household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point.

Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner.

15. Technical specifications

Model:	BF 600
Dimensions:	30 x 30 x 2.1 cm
Weight:	1.6 kg
Measurement range:	5-180 kg
Scale interval d:	The display can be read in increments of 0.1 kg.
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 0.4 kg (several measurements in direct succession on the same scale with position of scale and person the same as far as possible).
Absolute precision:	In comparison with a calibrated weight, the measured value is +/-1% +0.1 kg. E.g. at 40 kg this corresponds to +/- 0.5 kg; at 100 kg this corresponds to +/-1.1 kg.
Data transfer via Bluetooth ® low energy technology	The product uses Bluetooth ® low energy technology, frequency band 2400 - 2483.5 MHz, transmission power max. -0.93 dBm.

<p>A list of compatible smartphones, information on software and details on the devices can be found at the following link: connect.beurer.com</p>	
<p>System requirements for the "beurer HealthManager Pro" / "beurer HealthManager" app</p>	<p>iOS ≥ 12.0, Android™ ≥ 8.0 Bluetooth® ≥ 4.0.</p>

Subject to technical changes.

We hereby confirm that this product complies with the European RED Directive 2014/53/EU. The CE Declaration of Conformity for this product can be found under:

www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php.

16. Warranty/service

Further information on the warranty and warranty conditions can be found in the warranty leaflet supplied.

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The *Bluetooth*[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.

Beurer GmbH • Söflinger Str. 218 • 89077 Ulm (Germany) • www.beurer.com
www.beurer-gesundheitsratgeber.com • www.beurer-healthguide.com