beurer

BF 105 Body Complete

DE Diagnosewaage Gebrauchsanweisung	2	IT Báscula de diagnóstico Istruzioni per l'uso	69
EN Diagnostic bathroom scale Instructions for use	19	TR Diyagnoz terazisi Kullanım kılavuzu	86
FR Pèse-personne impédancemètre Mode d'emploi	36	RU Диагностические весы Инструкция по применению	102
ES Báscula de diagnóstico Instrucciones de uso	53	PL Waga diagnostyczna Instrukcja obsługi	120

CE

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Read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Contents

1. Getting to know your device	19	10. Evaluating the results	29
2. Signs and symbols	20	11. Transferring measurements	31
3. Warnings and safety notes	21	12. More information	31
4. Information	22	13. Cleaning and maintaining the device	32
5. Unit description	23	14. Disposal	32
6. Initial use	23	15. What if there are problems?	32
7. Set-up with app	24	16. Technical specifications	33
8. Set-up without app	25	17. Warranty / service	34
9. Taking measurement	26		

Included in delivery

- Diagnostic scale
- Removable operating unit with display
- Wall bracket
- Mounting kit for wall fastening
- 3 x 1.5V type AA (scale)
- 3 x 1.5V type AA (operating unit)
- USB cable
- These instructions for use and *Bluetooth*[®]/USB quick guide

Dear customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty, baby and air.

Please read these instructions for use carefully and keep them for later use. Be sure to make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer team

1. Getting to know your device

Functions of the device

The purpose of this digital scale is weighing and calculating your personal fitness data. It is intended for private use.

The scale is equipped with the following functions, which can be used by up to 10 different people:

- Body weight measurement
- Calculating body fat percentage,
- Body water content,
- Muscle percentage,
- Bone mass and
- Basal metabolic rate and active metabolic rate.

The scale also includes the following additional features:

- Innovative menu navigation on large display (in 5 languages),
- The display can be left on the scale, positioned separately or mounted on the wall. A stand kit is optionally available,
- Switch between kilograms “kg”, pounds “lb” and stones “st”,
- Automatic switch-off function,

- Low battery indicator,
- Automatic user recognition,
- Storage of the last 30 measurements for 10 users if the data cannot be transferred to the app.

2. Signs and symbols

The following symbols appear in these instructions for use:

	WARNING Warning notice indicating a risk of injury or damage to health.		
	IMPORTANT Safety note indicating possible damage to the device/accessory.		
	Note Note on important information		
	Observe the instructions for use		
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE		
	Dispose of packaging in an environmentally friendly manner		
	This product satisfies the requirements of the applicable European and national directives.		
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.		
	Manufacturer		
	The products demonstrably meet the requirements of the Technical Regulations of the EAEU.		
	Do not load the scale beyond 180 kg/396 lb/28 st.		Danger of slipping: do not step on the scale with wet feet.
	Danger of tipping: position yourself in the centre of the weighing surface.		Place the scale on an even surface. No carpet.

3. Warnings and safety notes

WARNING

- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- Do not use during pregnancy.
- Never step onto the outer edge of the scale with one foot as this could cause the scale to tip!
- Keep packaging material away from children (risk of suffocation).
- Warning, do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- **⚠ Choking hazard!** Small children could swallow and choke on batteries. Store batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- **⚠ Risk of explosion!** Do not throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device will not be used for some time, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

General notes

- The device is suitable for private use only, and is not intended for medical or commercial purposes.
- Please note that measuring tolerances are possible for technical reasons as this scale is not calibrated for use in a professional medical context.
- In cases of differing measurements (between the scale and app), only use the measured values displayed on the scale.
- The scale's maximum capacity is 180 kg (396 lb/28 st). The results for weight and bone mass are displayed in 100 g increments (0.2 lb).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- Calorie requirement is indicated in steps of 1 Kcal.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm".
- Place the scale on an even, hard surface; a hard surface is crucial to accurate measurements.
- Repairs may only be carried out by Beurer Customer Services or an authorised supplier. Before submitting a complaint, please check the batteries first and replace them if necessary.

Storage and maintenance

The accuracy of the measurements and service life of the device depend on its careful handling:

IMPORTANT

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when not in use.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Do not press the buttons violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile telephones).

4. Information

The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content in a matter seconds, using a current that cannot be felt, is completely harmless and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

Moreover, this scale uses two different frequencies to be able to determine impedance even better. These different frequencies affect the cell membrane and body water differently. The electrical contact to the body is also enhanced using 4 hand electrodes and 4 foot electrodes to achieve a higher measurement accuracy and better analysis. This sensitive analysis technology means readings may change slightly when measurements are repeated. This is caused by changes during the impedance measurement (e.g. different electrical contact to the electrodes by the hands and feet, different distribution of water in the body).

Carry out the analysis as follows:

- Upper body: Hold the control unit in your hands and step onto the scale wearing your shoes.
- Lower body: Step onto the scale barefoot. Do not hold the control unit this time (e.g. leave it on the scale).
- Upper body and lower body: Hold the control unit in your hands and step onto the scale barefoot.

If you are holding the control unit in your hands when performing the analysis, the weight of the unit is automatically deducted.

Please note that the values calculated by the diagnostic scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

General tips

- As far as possible, weigh yourself at the same time every day (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be made when barefoot; the soles of the feet may be lightly dampened where appropriate.
- Results may be unsatisfactory if the soles of the feet are completely dry or have a large amount of hard skin since conductivity will be impaired.
- Stand up straight and still whilst weighing yourself.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in your body can disperse.
- Remember that only the long-term trend is important. Short-term changes in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.


Limitations

When calculating body fat and the other values, unusual and implausible results may occur in the case of:

- Children under 10 years of age,
- Performance athletes and body-builders,
- Persons with a fever, persons undergoing dialysis, persons with oedema symptoms and persons suffering from osteoporosis,
- Persons taking cardiovascular medicine (affecting the heart and vascular system),
- Persons taking vasodilative or vasoconstrictive medication,
- Persons with significant anatomical abnormalities in their legs compared with their overall body height (leg length significantly shorter or longer).

5. Unit description

Overview



	1. Foot electrodes, left foot
	2. Control unit with display
	3. Operating status display 
	4. Foot electrodes, right foot

Control unit with display

	5. Contact electrodes, left hand
	6. Contact electrodes, right hand
	7. set Set button, confirm the entry
	8. ▲ Menu: up button
	9. ▼ Menu: down button
	10.  ON button
	11. ▲ Setting: up button
	12. ▼ Setting: down button

6. Initial use

Inserting the batteries


Insert the batteries into the scale and the control unit the correct way round:	Control unit	Scale
		

Positioning the scale

Place the scale on a flat and stable surface. A solid surface is a prerequisite for a precise measurement.


Setting up the control unit

You have three options:

<p>a) Attach to scale</p> 	<p>b) Set down on a flat surface</p> 
--	---

c) Wall fastening

- Use the supplied mounting kit for mounting and mount the wall bracket as shown. Use a suitable wall plug for the surface if necessary.
- Use the wall bracket as a drilling template to mark the holes.

You can also use the double-sided adhesive strip to attach the wall bracket to the wall. The wall surface must be clean and dry. Suitable for affixing to tiles, glass, plastic and many other surfaces.


- Take the backing off the adhesive strip on the back of the wall bracket.
- Push the wall bracket firmly against the wall.
- Check the bracket is securely attached.

Wait for at least one hour after mounting the wall bracket using the adhesive strip before fitting the control unit on the bracket!

To make fitting easier, first hold the control unit above the wall bracket against the wall and then slide the control unit down.

You can swivel the control unit on the wall bracket to adjust the angle. Locking the control unit ensures that it cannot be removed from the wall bracket when swivelled.

Swivel the control unit all the way to the front before taking it off the wall bracket!

As the perfect complement for even more comfort and as a highlight in the fitness/wellness range, we recommend the Professional Set consisting of a compatible aluminium column + glass base.

Please visit www.beurer.com or contact the corresponding service address (according to the service address list) in your country.

Designation	Item number and/or order number
Profi-Set	748.11

7. Set-up with app

To calculate your body fat percentage and additional body values, you must save the personal user data on the scale.

The scale has 10 user memory spaces which allows, for example, you and your family members to save personal settings.

It is also possible to activate users via other mobile end devices on which the “beurer HealthManager Pro“ / “beurer HealthManager“ App has been installed or by changing the user profile in the app (see the app settings).

Keep the smartphone close to the scale in order to maintain an active *Bluetooth*® connection during initial use.

- Activate *Bluetooth*® in the smartphone settings.
- Insert the batteries and place the scale on a solid surface (“CONNECT WITH APP” will flash on the scale’s LCD).
- Download the free “beurer HealthManager Pro“ / “beurer HealthManager“ App from the Apple App Store (iOS) or from Google Play (Android™).
- Start the app and follow the instructions.
- Enter the BF 105 related data.

The following settings must be set or entered in the “beurer HealthManager Pro” / “beurer HealthManager” App:

User data	Set values
Initials	Max. 3-digit (A-Z, 0-9)
Gender	Male (♂), female (♀)
Height	100 to 250 cm (3' 3.5" to 8' 02")
Date of birth	Year, month, day
Activity level	1 to 5
Target weight	kg/lb/st (for progress graph)

Activity levels

The medium-term and long-term perspective is decisive when selecting the activity level.

Activity level	Physical activity
1	None.
2	Low: little and light physical exertion (e.g. going for walks, light gardening work, gymnastic exercises).
3	Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.
4	High: physical exertion, at least 4 to 6 times a week for 30 minutes each.
5	Very high: intensive physical exertion, intensive training or hard physical work, daily for at least 1 hour.

The app assigns the next free user memory space (person 1 – person 10) to you.

- Carrying out the measurement: For automatic personal recognition, the first measurement must be assigned to your personal user data. The app will prompt you to do this. See section „9. Taking measurement“ for how to do this.

8. Set-up without app

We recommend using the app for the initial use and implementing all settings. See section „7. Set-up with app“. However, this can also be carried out on the device itself.

- Insert the batteries and place the scale on a solid surface (“CONNECT WITH APP” will flash on the scale’s LCD).
- Specify device/time/user settings on the control unit. To do so, use the [▽] and [△] buttons and confirm with [SET].

a. Specifying device settings

[SET] → “SET_UNIT ”

Sequence	Set values
Language	DE EN FR ES IT
Measurement unit	kg – lb – st
Time	24 h – 12 h
Guest mode	ON_OFF
Bluetooth®	ON_OFF
	 When set to “OFF”, a Bluetooth® connection is only possible if the handheld unit is switched on.

b. Setting the time

[SET] → “SET_TIME”

Sequence	Set values
Date	Year – Month – Day
Time	Hours – Minutes

c. Setting user data

[SET] → "SET_USER"

User data	Set values
Initials	Max. 3-digit (A-Z, 0-9)
Gender	Male (♂), female (♀)
Height	100 to 250 cm (3' -03" to 8' -02")
Date of birth	Year, month, day
Activity level	1 to 5 (see section 7)
Target weight	kg/lb/st, in kg (for progress graph)
User PIN	Assigned automatically for <i>Bluetooth</i> ® use at a later time

– Carrying out the measurement: For automatic personal recognition, the first measurement must be assigned to your personal user data. See section „9. Taking measurement“ for how to do this.

Logging in to the app at a later time


If you have carried out the initial use of the scale without using the app, you can set up your user on the app as follows.

- Download the free “beurer HealthManager Pro“ / “beurer HealthManager“ App from the Apple App Store (iOS) or from Google Play (Android™).
- Start the app and follow the instructions.

9. Taking measurement

Switch on scale/measure weight

– Press the  button on the control unit.
0.0 kg appears in the display.


– The scale is then ready to take a weight measurement. To measure your weight, step onto the scale with your weight distributed evenly on both feet.

 For a new measurement, the scale assigns the measurement to the user whose most recently saved measurement falls within +/3kg of this. If automatic user assignment is not successful, only the weight without initials is displayed. This happens if no definite assignment is possible or the measured value is outside +/3kg. In this case, you should select the user manually.

Manual user selection:

Select your user by pressing the right  button. If necessary, press the button several times.

Press the [SET] button or wait three seconds.

 If you activate the user “GUEST”, you will need to re-enter the user data before every measurement. In “GUEST” mode, the user data and measured values are not stored. You have the option to switch “GUEST” mode on or off, see section 8a.

Diagnosis

With the control unit switched on, stand on the scale barefoot holding the control unit in your hands. Hold the control unit with the underside of your thumbs on the front electrodes and with your fingers pressed firmly against the back electrodes. When you are holding the control unit in your hands, the weight of the unit is automatically deducted.


Stand still with your weight distributed equally on both feet. Make sure there is no skin contact between your two feet, legs, calves and thighs.

Following the weight measurement, impedance will be measured for the analysis: the upper body using the 4 electrodes on the control unit, and the lower body using the 4 electrodes on the scale. This is indicated by the symbol showing current passing through the body.

The overview is displayed.

The measured value is automatically saved when switching off.

Note: The measurement will not be saved if the following applies:

The user is not automatically detected. This happens if your measured weight differs from the previous measurement by +/3 kg or if another user has similar values (+/3 kg) to yours. In this case, press the [△] button to select the user manually.

The user is standing on the scale with shoes on and is not holding the handheld unit.

Analysis of measured values

The results will be analysed automatically after taking the measurement or after selecting the user manually on the control unit.

Graphic representation


a) Weight

– Press the menu button [▼] on the control unit to switch to the individual sub-menus or the memory.

The trend for your body weight is shown from measurement to measurement. Later measurements are added on the right. The reference line shows your target weight.


You can change your target weight in the user data at any time.

b) Body fat

– Use the [▽, △] buttons to switch to showing body fat values.

The trend for your body fat is shown from measurement to measurement. The reference line represents the average value for all body fat measurements; the bars show whether the values have gone up or down.


Upper body, lower body

a) Body fat


- Press the [▼] button on the control unit to switch to the detailed view for the upper/lower body.

The current body fat measurement is displayed for the whole body and is subdivided into upper body/lower body.

b) Muscle

- Use the [▽, △] buttons to show the detailed muscle mass view.
-  Only the current measurements can be shown, not the old stored measured values.


Mem Show stored measurements

- Press the [▼] button on the control unit to switch to the **memory**.

The overview is displayed. The measurement date and the memory space number are shown at the top, e.g. (03-03) means the third space = the last of a total of 3 occupied memory spaces.

- Use the [▽, △] buttons to cycle through the memory spaces.


Average values

- Press the [▼] button on the control unit to switch to the **average**.


All average values in the stored measurements are displayed.

-  The dates of the first and the last measurement are shown at the top.

Kcal Calorie consumption

- Press the [▼] button on the control unit to show **Kcal**. The theoretically calculated basal metabolic rate BMR and the active metabolic rate AMR are shown in Kcal.

Switching off the scale (standby mode)

- Wait until the display switches itself off or press the [U] button on the control unit to switch the display off manually.

-  The date and time remain visible.

10. Evaluating the results

Body fat percentage

Body fat values are shown in the display and their analysis is shown in the bar chart. The subsequent body fat percentages serve as a guideline (for further information ask your doctor).

Male

Age	Low	Normal	High	Very high
10-14	<11 %	11-16 %	16.1-21 %	>21 %
15-19	<12 %	12- 17 %	17.1-22 %	>22 %
20-29	<13 %	13-18 %	18.1-23 %	>23 %
30-39	<14 %	14-19 %	19.1-24 %	>24 %
40-49	<15 %	15-20 %	20.1-25 %	>25 %
50-59	<16 %	16-21 %	21.1-26 %	>26 %
60-69	<17 %	17-22 %	22.1-27 %	>27 %
70-100	<18 %	18-23 %	23.1-28 %	>28 %

Female

Age	Low	Normal	High	Very high
10-14	<16 %	16-21 %	21.1-26 %	>26 %
15-19	<17 %	17-22 %	22.1-27 %	>27 %
20-29	<18 %	18-23 %	23.1-28 %	>28 %
30-39	<19 %	19-24 %	24.1-29 %	>29 %
40-49	<20 %	20-25 %	25.1-30 %	>30 %
50-59	<21 %	21-26 %	26.1-31 %	>31 %
60-69	<22 %	22-27 %	27.1-32 %	>32 %
70-100	<23 %	23-28 %	28.1-33 %	>33 %

Values will often be lower for fitness enthusiasts. Depending on the sporting activities performed, the training intensity and the physical constitution, the results may still be below the stated standard values. Please note, however, that there may be a risk to health in the case of extremely low values.

Body water content

The body water content in % is normally within the following ranges:

Male

Age	Poor	Good	Very good
10-100	<50 %	50-65 %	>65 %

Female

Age	Poor	Good	Very good
10-100	<45 %	45-60 %	>60 %

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle percentage

The muscle percentage is normally within the following ranges:

Male

Age	Low	Normal	High
10-14	<44%	44–57 %	>57 %
15-19	<43%	43–56 %	>56 %
20-29	<42%	42–54 %	>54 %
30-39	<41%	41–52 %	>52 %
40-49	<40%	40–50 %	>50 %
50-59	<39%	39–48 %	>48 %
60-69	<38%	38–47 %	>47 %
70-100	<37%	37–46 %	>46 %

Female

Age	Low	Normal	High
10-14	<36%	36–43 %	>43 %
15-19	<35%	35–41 %	>41 %
20-29	<34%	34–39 %	>39 %
30-39	<33%	33–38 %	>38 %
40-49	<31%	31–36 %	>36 %
50-59	<29%	29–34 %	>34 %
60-69	<28%	28–33 %	>33 %
70-100	<27%	27–32 %	>32 %

Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. Bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Please note that this scale does not identify the total bone mass, but rather just the mineral content of the bone (without water content and without organic substances). It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

Important

Please do not confuse bone mass with bone density. The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. lying in bed 24 hours a day). This value is primarily dependent on weight, height and age. It is displayed on the diagnostic scale in the unit kcal/day and calculated using the scientifically recognised Harris-Benedict equation.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic scale calculates this by means of the activity level (1–5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink. If less energy is taken in than is used over a long period of time, the body takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Temporal context of the results

Please note that only the long-term trend is significant. Brief deviations in weight within a few days are normally caused by loss of fluids.

The interpretation of the results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as on the length of time over which these changes occur. Rapid changes within the scope of a few days are to be considered separately from medium-term changes (in the scope of weeks) and long-term changes (months).

It can be said as a basic rule that short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may also relate to the fat percentage and the muscle percentage.

- If your weight drops in the short term but your body fat percentage increases or stays the same, you have only lost water, for example following a training session, visit to the sauna or a diet aimed only at fast weight loss.
- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have built up valuable muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working – you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water and muscle percentages should not be totalled (muscle tissue also contains components made of body water).

11. Transferring measurements

a) Using the app

If the app is not open, the newly assigned measurements are saved on the scale. A total of 30 measurements per user can be saved on the scale. The saved measurements are transferred automatically to the app when you open the app within the *Bluetooth®* range. It is possible to automatically transfer data when the scale is switched off.

b) By USB

- You will find the Beurer software download and installer in the download area at <http://www.beurer.com>
- Follow the instructions for software installation.
- Establish USB connection.

12. More information

Deleting scale data

a) Delete user with all settings and saved values

- After switching on select your user by pressing the [\triangle] button. If necessary, press the button several times and confirm with [SET].
- When  appears in the display, hold down the [SET] button for **10** seconds to delete the user.
- Use the [∇ , \triangle] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the user, confirm “NO” using the [SET] button.

All saved values and settings for this user will be deleted.

b) Delete individual saved value for a user

- Select a measured value in the measurements memory “Mem” and briefly press the [SET] button on the control unit to delete the stored value displayed.
- You can also do this immediately after the measurement is taken.
- Use the [∇ , \triangle] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete the value, confirm “NO” using the [SET] button.

c) Delete all stored values for a user

- Open the measurements memory “Mem” and press the [SET] button on the control unit for around 5 seconds to delete all values.
- Use the [∇ , \triangle] buttons to select “YES” and confirm the prompt using the [SET] button. If you do not want to delete, confirm “NO” using the [SET] button.

Changing the batteries

The control unit has a low battery indicator. The  symbol appears in the display when the control unit batteries are flat, and the  symbol appears if the scale batteries are flat. In this case, the batteries should be replaced. Connect the scale to the app again. The system time is updated.

**NOTE:**

- When changing the batteries, use batteries of the same type, make and capacity.
- Do not use rechargeable batteries.
- Use batteries free from heavy metals.

13. Cleaning and maintaining the device

The device should be cleaned from time to time.

This should be done using a damp cloth and, if necessary, a small amount of detergent.

**IMPORTANT**

- Never use abrasive solvents or cleaning products!
- Never submerge the device in water!
- Do not clean the device in a dishwasher.

14. Disposal


Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of batteries.

The codes below are printed on batteries containing harmful substances.

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

For environmental reasons, do not dispose of the device in the household waste at the end of its service life.

Dispose of the device at a suitable local collection or recycling point in your country.

Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner.

**15. What if there are problems?**


If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
ERROR:_OVERLOAD	Maximum weight capacity of 180 kg exceeded.	Loads must be limited to 180 kg.
LOW_BATTERY	Batteries are nearly flat.	Replace the batteries (see section 12).
MEMORY_FULL	30 measurements have been stored.	Transfer the measured values.
ERROR:_CONTACT	The electrical resistance between the electrodes and the soles of the feet is too high (e.g. with very hard skin).	Repeat the measurement barefoot. Dampen the soles of your feet slightly. If necessary, remove hard skin on the soles.
ERROR:_DATA	The BIA was successful, but the fat percentage is outside the measurable range (less than 5% or greater than 50%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary. If necessary, change your set activity level.
Empty display	The batteries in the control unit are completely flat.	Replace the batteries (see section 12).
	The batteries have not been inserted in the control unit correctly.	Check the polarity is correct, (see page 23).
Only weight is displayed, no diagnosis	Unknown measurement as it is outside the user assignment limit or a unique assignment is not possible.	Select user on the scale using the “ Δ ” button or repeat user assignment in the app. See section 9.

No measurement possible

Possible faults	Solution
The control unit was not activated before stepping onto the scale.	Switch the control unit on before taking the measurement. If necessary, wait 30 seconds and then repeat the measurement.
Radio connection between scale and control unit interrupted.	Re-insert the batteries in the scale and control unit. Ensure a sufficient distance from strong electromagnetic fields.
The control unit is too far from the scale.	Reduce the distance between the scale and the control unit to a maximum of 2 metres.
The batteries in the scale are flat or have been inserted incorrectly.	Change the batteries and check the polarity is correct.

Incorrect weight measurement

Possible faults	Solution
Incorrect zero point of scale.	Just before taking the measurement, activate the scale with the foot until the operating status  appears in the display. Wait approx. 5 seconds and repeat the measurement. For your information: because you step on the scale straight away in normal use, the scale saves the zero point when the scale is switched off (when the operating status indicator  goes out).
Scale is on a carpet.	Place the scale on a flat, solid surface. Stand on the outside of the scale as far as possible without it tipping.
Incorrect weight is displayed.	Stand as still as possible.
Incorrect weight is displayed.	Wait until the scale switches itself off. Then repeat the measurement.

No *Bluetooth*[®] connection (symbol is missing)

Possible faults	Solution
Radio connection	The range in the open is approx. 10 m. Walls and ceilings reduce the range. Other radio waves may interfere with the transmission. Therefore, do not position the scale near devices such as WIFI routers, microwaves, induction hobs etc. and, ideally, switch off any other <i>Bluetooth</i> [®] devices that are within range.
<i>Bluetooth</i> [®] connection	Close the app completely (including in the background). Switch <i>Bluetooth</i> [®] off and back on again. Switch the smartphone off and back on again. Remove the batteries in the scale for a short period and reinsert them. See FAQs at www.beurer.com .
Connection set-up taking too long	Faster <i>Bluetooth</i> [®] interval when the scale is switched on.

16. Technical specifications

Scale dimensions	330 x 330 x 41 mm
Scale batteries	3 x 1.5V type AA
Control unit batteries	3 x 1.5V type AA
Wireless transmission (scale to control unit)	433.92MHz; Maximum transmission power (ERP): -21.5dBm
Measurement range	5 kg to 180 kg

Scale interval d:	The display can be read in increments of 0.1 kg.
Repetition accuracy:	The measuring tolerance for repeated measurements is +/- 0.4 kg (several measurements in direct succession on the same scale with position of scale and person the same as far as possible).
Absolute precision:	In comparison with a calibrated weight, the measured value is +/-1% +0.1 kg. E.g. at 40 kg this corresponds to +/- 0.5 kg; at 100 kg this corresponds to +/-1.1 kg.
Data transfer via <i>Bluetooth</i> [®] wireless technology	This product uses <i>Bluetooth</i> [®] low energy technology, frequency band 2.402 – 2.480 GHz, transmission power max. 0 dBm, compatible with <i>Bluetooth</i> [®] 4.0 smart-phones/tablets
A list of compatible smartphones, plus information and software for the devices can be found at the following link: connect.beurer.com	
System requirements for the “beurer HealthManager Pro” USB uploader and for the “beurer HealthManager” PC software	Windows 8.1 and above USB 2.0 (type A) and above
System requirements for the “beurer HealthManager Pro” / “beurer HealthManager” App	iOS ≥ 12.0, Android™ ≥ 8.0 <i>Bluetooth</i> [®] ≥ 4.0

Subject to technical changes.

We hereby confirm that this product complies with the European RED Directive 2014/53/EU. The CE Declaration of Conformity for this product can be found under:
www.beurer.com/web/we-landingpages/de/cedeclarationofconformity.php.

17. Warranty / service

Beurer GmbH, Söflinger Straße 218, 89077 Ulm, Germany (hereinafter referred to as “Beurer”) provides a warranty for this product, subject to the requirements below and to the extent described as follows.

The warranty conditions below shall not affect the seller’s statutory warranty obligations which ensue from the sales agreement with the buyer.

The warranty shall apply without prejudice to any mandatory statutory provisions on liability.

Beurer guarantees the perfect functionality and completeness of this product.

The worldwide warranty period is 5 years, commencing from the purchase of the new, unused product from the seller.

The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.

German law shall apply.

During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, Beurer shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.

If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached “International Service” list of service addresses.

The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.

A warranty claim shall only be considered if the buyer can provide Beurer, or an authorised Beurer partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this warranty:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by Beurer;
- damage that arises during transport between manufacturer and customer, or between service centre and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the warranty period under any circumstances.

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The *Bluetooth*[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android[™] is a trademark of Google LLC.