




FT 100 Non-contact clinical thermometer


 german
 engineering


Measures without skin contact - hygienic and safe



Distance sensor (LED/acoustic signal)



LED temperature alarm (green/yellow/red)



Measurement of forehead, object and room
 temperature in seconds

60 memory spaces (forehead measuring)

°C/°F can be set

Display illumination (On/Off)

Sound control (On/Off)

Automatic switch off

Incl. storage box

With battery changing indicator

Incl. 2 x 1.5V batteries AAA

Product measurements (LxWxH):

ca. 199 x 68.2 x 40.9 mm

Product weight: approx. 82 g

Medical device

5 year guarantee

Sales unit: 2 / Shipping carton: 8 x 2

EAN no.: 4211125 79505 4

Item no.: 795.05



< 37,5°C
 everything is OK



37,5-37,9°C
 elevated temperature



≥ 38°C
 fever alarm