



PO 13
Kids

beurer
babycare

PO 13 Kids Pulse Oximeter for Kids



Determines the oxygen saturation (SpO₂) and heart rate (pulse)



Very simple and painless measurement within 8 seconds

Displays measurements within seconds without testing the child's patience

Adapted to suit children's fingers*

Small and light for use at home and on the move

Particularly suitable for children with
- cardiovascular diseases
- respiratory diseases

Easy-to-read colour display

Colour display with 4 viewing perspectives

Adjustable display brightness

The display is oriented automatically (portrait, landscape)

Graphic pulse display

Low battery indicator

Automatic switch-off

Incl. strap, belt pouch and screw driver

Medical device according to new Medical Device Regulation (EU) 2017/745 (MDR)

Incl. 2 x 1.5 V AAA batteries

Product dimensions: 59 x 37 x 35 mm

Product weight: approx. 27 g (without batteries)

5-year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 421112560095 2

Item no.: 600.95



* Recommended finger width: 0.7 – 1.1 cm,
Recommended weight: 10 – 40 kg