



BC 51

beurer
medical

BC 51 – Wrist blood pressure monitor



positioning
indicator

Accurate measurement results due to correct measurement position at heart level



inflation
technology

Fast and convenient measurement, even during pumping up



Fully automatic blood pressure and pulse measurement on the wrist

With cuff for wrist circumferences from 12.5 - 21.0 cm

Large, perfectly readable XL display

Risk indicator: Classification of measurements with coloured scale

Arrhythmia detection: Warning in case of possible heart rhythm disturbances

Extra slim design – the ideal companion, even on the go

Two user memories for 120 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 2 x 1,5V AAA batteries and practical storage case

Medical device

Product measurements (LxWxH): 95 x 70 x 80 mm

Product weight: approx. 105 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65060 5

Item no.: 650.60