



## IH 51 Nebuliser



ultraquiet



Relaxed usage thanks to whisper-quiet operation

home and away



Small and light - ideal for at home and on the road

selfcleaning function



Hygienic due to self-cleaning function



Innovative nebulisation using vibrating membrane technology

For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory diseases, etc.

- for fast and effective inhalation

Particle size (MMAD): 3,89  $\mu\text{m}$

High nebulisation capacity (Low nebulisation rate: approx. 0.2 ml/min, High nebulisation rate: approx. 0.4 ml/min) - short inhalation time

Accessories: Mesh atomiser with medication container, mouthpiece, adult mask, child mask, storage bag and USB-C cable

Medical device according to new Medical Device Regulation (EU) 2017/745 (MDR)

Product measurements (LxBxH): 46 x 62 x 112 mm

Product weight: approx. 0,093 kg

Yearpack: 110.124

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN No.: 4211125 103859

Item no.: 103.85