




IH 16 Nebuliser



Space-saving - small and compact design



Easy & hygienic cleaning thanks to removable elements



Holistic treatment of the airways for adults and children

A proven way to treat acute and chronic diseases, for respiratory care and to prevent colds

Compressed-air technology produces the finest particles, which transport the active ingredients deep into the lungs

Particle size (MMAD): 2,07 µm

Medication flow rate: approx. 0.3 ml/min

Atomiser set: Nebuliser bottle, mouthpiece, nosepiece, adult mask, child mask, compressed air tube, 5 filter

Medical device according to new Medical Device Regulation (EU) 2017/745 (MDR)

Operating pressure: approx. 0.35 - 0.8 bar

Mains operation: 230 V~ / 50 Hz / 230 VA

Product measurements (LxBxH): 138.6 x 116 x 87.8 mm

Product weight: approx. 0.85 kg

Yearpack: 110.125

5 year guarantee

PZN: 19362160

HMV: in application

Sales unit: 4 / Shipping carton: 4 x 2

EAN No.: 4211125 103514

Item no.: 103.51

