



## BM 23 Upperarm Blood pressure monitor



Easy to use via one single button



With USB-C connection for optional mains operation



Scan & Save – the simple way  
to save blood pressure  
values in the app



Fully automatic blood pressure and pulse measurement on the upper arm

Deflation Technology

Cuff position check: For ensuring the cuff is correctly positioned on the upper arm

Arrhythmia detection: detects and warns of possible cardiac arrhythmias

Risk Indicator: Classification of measurements with coloured scale

Automatic switch-off

Battery changing indicator

With cuff for upper arm circumferences from 22 - 42 cm

Incl. 4 x 1,5V AAA batteries

Medical device according to new Medical Device Regulation (EU) 2017/745 (MDR)

Product measurements (LxWxH): 120 x 90 x 40 mm

Product weight: approx. 318 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 6 x 4

EAN no.: 4211125 60493 6

Item no.: 604.93