

EM 89  
Heat


beurer  
medical



## EM 89 Digital TENS/EMS unit with heat function

german  
engineering

PAIN  
free



TENS  
technology

Transcutaneous electrical  
nerve stimulation (TENS)




EMS  
technology

Electrical muscle stimulation (EMS)



heat  
function

TENS/EMS and heat function  
can be used in parallel and/or  
individually



4-in-1 Function: TENS, EMS, Massage, Heat

4 separately adjustable channels with 8 self-adhesive electrodes

64 programs within the categories TENS, EMS and Massage

6 free programmable programs (TENS/EMS): Frequency, pulse width, on time, off time

"Favourite programme": Supports easy and quick access to the favoured programme

Backlight illuminated, easy readable display

Safety switch-off

Output voltage: max. 100 Vpp (500 Ohm load)

Output current: max. 200 mA<sub>pp</sub> (500 Ohm load)

Pulse duration: 50-450 µs

Output frequency: 1-150 Hz

Pulse form: dual-phase square pulses

With detailed instruction manual (inclusive application possibilities)

Lithium-ion rechargeable battery (4000 mAh, 3.7 V)

Convenient storage pouch included

Medical product

Size: ca. 142 x 159 x 53 mm

Weight: ca. 341g

5 year guarantee

Sales unit: 2 / Shipping carton 4 x 2

EAN no.: 4211125102777

Item no.: 102.77

102.77\_EM89\_2023-11-09\_02\_DS\_DS\_EN\_BEU Subject to errors and changes



Drug Free  
Pain Relief