



BM 58

beurer  
medical

## BM 58 Blood pressure monitor



Easy to use via sensor touch buttons



Clear, easy-to-read XL black display



For innovative networking between the blood pressure monitor and "beurer HealthManager" PC software



Fully automatic blood pressure and pulse measurement on the upper arm

Arrhythmia detection:

Warning in case of possible heart rhythm disturbances

Two user memories for 60 measured values each

Average value of all stored measurements and morning and evening blood pressure for the last 7 days

With cuff for upper arm circumferences from 22 - 30 cm

Cuff holder

Can be used with mains part

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

USB cable, 4 x 1.5V AA batteries and storage pouch included

Medical device

Mains part available separately: Item no. 071.95

Product measurements (LxWxH): 100 x 150 x 58 mm

Product weight: approx. 508 g (with cuff, without batteries)

XL-cuff (30 - 42 cm) available separately: Art.-No. 163.246

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65516 7

Item no.: 655.16



**Sensor touch**

- Easy to use via sensor touch buttons




**Mains part ready**

- For optional operation of the device with separately available mains part



## CONNECTION



Innovative networking between the blood pressure monitor and the "beurer HealthManager" PC software with USB cable included - keep track of your measurements.



\* tested online application "beurer HealthManager" Cloud  
Can only be used with its full functionality with the „beurer HealthManager“ software.