



BM 27

beurer  
medical

## BM 27 Blood pressure monitor



universal  
cuff

With cuff for upper arm circumferences  
from 22-42 cm



cuff position  
control

Cuff position control\*: For ensuring  
the cuff is correctly positioned  
on the upper arm



Fully automatic blood pressure and pulse measurement on  
the upper arm

Clear, easy-to-read display

Risk indicator: classification of measurements with coloured  
scale

Arrhythmia detection: Warning in case of possible heart  
rhythm disturbance

Four user memories for 30 measured values each

Average of all measurements

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Battery changing indicator

Incl. 4x1,5 V batteries AA and storage bag

Medical device

Product measurements (LxWxH): 112 x 110 x 53 mm

Product weight: approx. 341 g (with cuff, without batteries)

5 year guarantee

Sales unit: 4 / Shipping carton: 6 x 4

EAN no.: 4211125 65818 2

Item no.: 658.18

### \*Cuff wrapping guide

