



EM 95 Bluetooth® - **EMS** homeSTUDIO Muscle stimulation device


 german
 engineering


Electric muscle stimulation (EMS)
 High End EMS device for home use



Cuffs for arms & legs



App with 20 predefined workouts and
 Virtual Coach



High End EMS device for home use

20 workouts in 2 categories:

Fitness & Power and Relax & Wellbeing

3.5" color TFT Display

4 channels

4 cuffs (2x arms and 2x legs)

40 programs

Handy belt clip

Li-Ion battery: 1170 mAh (4h run time)

USB charging cable including power adapter

Pouch

App with 20 predefined workouts and Virtual Coach incl. videos of
 all exercises


Compatible from iOS 9.0, Android™ 5.0, Bluetooth® 4.0

3 year warranty

VE: 2 / Transportkarton: 6x2

EAN no.: 4211125 66201 1

Item no.: 662.01



Practical trouser/belt clip for convenient use during the exercise

Four cuffs to stimulate the muscles in the upper arm and thighs


40 programs

Eight self-adhesive gel electrodes to stimulate the muscles (e.g. abdominal muscles)



Achieve your training targets more easily

Shape your body with the Beurer EM 95 Muscle stimulation device.



Just 3 steps to your success

- Choose and edit your workout
 - Positioning of electrodes and cuffs
 - Train with your Virtual Coach
-
- 20 predefined workouts and 50 exercise videos
 - Optimized for smartphone and tablet
 - With trainerfeedback and training progress in detailed statistics



CONNECTION

beurer
connect



Innovative networking between smartphone and muscle stimulation device - Select and create workouts at home or on the go in the app and transfer them to the device via Bluetooth®.



Syncs with

- iOS ≥ 9.0
- Android™ ≥ 5.0
- Bluetooth® ≥ 4.0



List of supported devices