



ETM TESTMAGAZIN - URTEIL	
Beurer BM 54	
SEHR GUT	94,1 %
Blutdruckmessgerät im Online-Einzeltest www.etm-testmagazin.de	
	04/2019



BM 54

beurer
medical

BM 54 Bluetooth® - Upper arm blood pressure monitor



Universal cuff: also suitable for larger upper arm circumferences from 22-44 cm



Risk indicator: classification of measurements with coloured scale



- Fully automatic blood pressure and pulse measurement on the upper arm
- Innovative networking between smartphone and blood pressure monitor via Bluetooth® low energy technology
- Arrhythmia detection: Warning in case of possible heart rhythm disturbances
- Large, easy-to-read XL display
- Two user memories for 60 measured values each
- Average value of all stored measurements and morning and evening blood pressure for the last 7 days
- Date and time
- Automatic switch-off
- Signal in case of application error
- Automatic pressure preselection and deflation
- Incl. 4 x 1,5V AAA batteries and practical storage pouch
- Medical device
- Product measurements (LxWxH): 139 x 94 x 48 mm
- Product weight: approx. 364 g (without batteries, with cuff)
- 5 year guarantee
- Sales unit: 4 / Shipping carton: 4 x 4
- EAN no.: 4211125 65512 9
- Item no.: 655.12



ETM TESTMAGAZIN - URTEIL	
Beurer BM 54	
SEHR GUT	94,1 %
Blutdruckmessgerät im Online-Einzeltest www.etm-testmagazin.de	
	04/2019



BM 54
beurer
medical



CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with

- iOS ≥ 10.0
- Android™ ≥ 5.0
- Bluetooth® ≥ 4.0



List of supported devices

* tested online application "beurer HealthManager" Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.