



EM 59 Digital TENS/EMS device with heat function



TENS
technology

Transcutaneous electrical
nerve stimulation (TENS)




EMS
technology

Electrical muscle stimulation (EMS)



heat
function

TENS/EMS and heat function
can be used in parallel and/or individually



PAIN
free

german
engineering

4-in-1 Function: TENS, EMS, Massage, Heat

2 separately adjustable channels with 4 self-adhesive
electrodes (electrodes included; size 45 x 45 mm)

64 pre-programmed applications (TENS/EMS/Massage)

6 individualizable programs (TENS/EMS)

"Doctor's-Function": ensuring ideal coordination to your
personal therapy programme

Timer function (adjustable between 5-100 minutes)

Inverse LCD with electrode position display for easy programme
selection and electrode placement

Ergonomic design, with belt clip

Safety switch-off

Output voltage: max. 100Vp-p (500 Ohm load)

Output current: max. 200 mA p-p (500 Ohm load)

Output frequency: 1-150 Hz

Pulse width: 50-450 µs per pulse

Pulse form: dual phase square pulses

With detailed instruction manual (inclusive application possibilities)

Rechargeable battery: 2000 mAh, 3,7V

Medical product

Product measurements (LxWxH): 139 x 66 x 26 mm (inkl. Clip)

Product weight: approx. 125 g (with belt clip, without batteries)

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4211125 64663 9

Item no.: 646.63

