

MG 55

beurer  
wellbeing



## MG 55 Klopf-Massagegerät

 german  
engineering



Deep-acting tapping massage



3 interchangeable massage attachments



For relaxing and loosening muscles

Continuously adjustable massage intensity

With on/off heat function

3 exchangeable massage attachments  
with rotary fixation for a safe support

Ergonomically shaped

Non-slip handle

Mains operation:

Input: 220-240 V~; 50/60 Hz; 20 W

Product measurements (LxWxH): 43.5 x 9.5 x 9.5 cm

Product weight: approx. 1000 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 64320 1

Item no.: 643.20