



FM 200
Achillomed®

beurer

FM 200 Achillomed® Achilles tendon massager

 german
engineering



To relieve problems with the Achilles tendon



Soothing massage of the
Achilles tendon that boosts circulation



6 rotating massage heads

2 massage speeds settings

With height and width adjustment

2 massage direction settings

Removable massage heads for easy cleaning

Illuminated LED operating display

Mains operation:

Input: 100-240 V~; 50/60 Hz; 0.45 A

Output: 15.0 V --- 1.0 A; 15.0 W

Product measurements (LxWxH): 26.5 x 33.5 x 13.4 cm

Product weight: approx. 2.5 kg

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64933 3

Item no.: 649.33



649.33_FM200_Achillomed_2020-10-23_02_DS_EN_BEU Subject to errors and changes

PhysioLine