

MG 180

beurer



# MG 180 MASSAGE GUN

## Deep tissue massager

 german engineering



4 Intensity levels



6 massage attachments for tailored use



Powerful and soothing trigger point massage



Targeted muscle massage for activation and regeneration

6 massage attachments for tailored use

4 intensity levels - up to 3.200 impulses per minute

With ergonomic handle for easy use

Can be used for all muscle areas

LED lights for displaying the battery state and operating setting

Ergonomically shaped

Powerful lithium ion battery for up to 3 hours of use

Product measurements (L x W x H): 17.5 x 7 x 24 cm

Product weight: approx. 1052 g

3 year guarantee

Sales unit: 1 / Shipping carton: 2

EAN no.: 4211125 64423 9

Item no.: 644.23



644.23\_MG180\_2021-04-09\_03\_DB\_EN\_BEU Subject to errors and changes

PhysioLine