




BF 105
Body Complete

beurer
wellbeing

BF 105 Body Complete Diagnostic bathroom scale



Complete body analysis: Accurate upper and lower body measurement with 8 electrodes



Innovative networking between smartphone and scale.

With free "beurer HealthManager Pro" /

"beurer HealthManager" app

Transfer using *Bluetooth*[®] low energy technology



Detailed body analysis with graphic evaluation



Determination of body weight, body fat, body water, muscle percentage and bone mass

Separate indication of upper and lower body measurement

for body fat and muscle percentage possible

Basic turnover and active turnover in kcal

Removable XXL display (suitable for wall mounting)

with blue LCD illumination

Time and date

Intuitive user guide in 5 languages (D/GB/F/ES/I)

Including USB cable for usage with PC

With brushed stainless steel electrodes

Change over to kg/lb/st

Capacity: 180 kg

Graduation:

- 0.1% for body fat, body water and muscle ratio

- 100 g for body weight and bone mass

5 activity levels

10 user memories with 30 measurement spaces each + guest mode

Overload indication

Button-On-Technology

Automatic switch-off

Incl. 6 x 1.5V batteries AA

Product measurements (LxWxH):

Scale: 33 x 33 x 4 cm Display: 21.2 x 7.7 x 3.8 cm

Product weight:

Scale: approx. 2.06 kg Display: approx. 260 g

5 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 74812 8

Item no.: 748.12

748.12_BF105_2021-09-21_01_DS_EN_BEU Subject to errors and changes





BF 105
Body Complete

beurer
wellbeing

BF 105 Body Complete Diagnostic bathroom scale



CONNECTION



Innovative networking between PC, smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via Bluetooth®



Syncs with
• iOS ≥ 12.0
• Android™ ≥ 8.0
• Bluetooth® 4.0



List of supported devices

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.