




FT 15/1 Express thermometer


 german
 engineering


Extra rapid measurement in approx. 10 seconds



Flexible measuring tip



Audible warning signal if temperature exceeds 37,8°C
 (fever alarm)



Digital EXPRESS thermometer

Audible signal when measuring stops

1 memory

Measurement accuracy $\pm 0.1^{\circ}\text{C}$ (between 35.5 and 42°C)

Measurement accuracy $\pm 0.2^{\circ}\text{C}$ (between 32 - 35.4°C and 42.1 - 42.9°C)

Extra large display

Waterproof

Can be disinfected

Without glass, mercury free

With storage case

Inkl. 1 x 1,5V Batterie LR41

Battery changing indicator

Automatic switch off

Indication of °C and °F

Delivered in sales display

Product measurements (LxWxH): 139 x 22 x 12 mm

Product weight: approx. 12 g (incl. Batterie)

Medical device

5 year guarantee

Sales unit: 20 / Shipping carton: 10 x 20

EAN no.: 4211125 79410 1

Item no.: 794.10