



# MG 185 MASSAGE GUN

## Deep tissue massager




9 intensity levels & 3 massage programmes



5 massage attachments for a tailored massage



Powerful and soothing trigger point massage



Powerful and soothing trigger-point massage for relaxing any tension

Strengthening and regeneration – quick and easy

5 massage attachments for tailored use

9 intensity levels & 3 massage programmes

- Intensity level 1: 1200 rpm or 20 Hz
- Intensity level 2: 1400 rpm or 23.3 Hz
- Intensity level 3: 1600 rpm or 26.6 Hz
- Intensity level 4: 1800 rpm or 30 Hz
- Intensity level 5: 2000 rpm or 33.3 Hz
- Intensity level 6: 2200 rpm or 36.6 Hz
- Intensity level 7: 2400 rpm or 40 Hz
- Intensity level 8: 2600 rpm or 43.3 Hz
- Intensity level 9: 3000 rpm or 50 Hz

Device is easy to use thanks to LED touch display

Very ergonomic with non-slip handle

Incl. case for practical storage

Powerful lithium ion battery for up to 5.5 hours of use

Product measurements (L x W x H): 15.5 x 6.5 x 26.0 cm (without attachments)

Product weight: approx. 1075 g (with attachments and adapter)

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 64533 5

Item no.: 645.33



645.33\_MG185\_2023-05-25\_05\_DB\_EN\_BEU Subject to errors and changes

PhysioLine