



IH 57 Nebuliser



vibrating
mesh

Nebulisation using vibrating membrane technology



intelligent
aerosol
technology

Intelligent aerosol technology –
effective application due to nebulization stop
during exhalation



lithium-ion
battery

Lithium ion battery with
a battery working time
of approx. 3.5 hours



For the treatment of the upper and lower respiratory tract

For use in colds, asthma, respiratory diseases, etc.
- for fast and effective inhalation

High proportion of respirable particles (MMAD): 4.52 µm

High nebulisation capacity (Low nebulisation rate: ≥0.2 ml/min
High nebulisation rate: ≥ 0.4 ml/min) - short inhalation time

Flexible rotation up to 45°

Ideal for children and suitable for travelling

Accessories:

Mesh atomiser with medication container, mouthpiece,
adult mask, child mask, storage bag and USB charging cable

Hygienic preparation due to self-cleaning function

Ultra-quiet (approx. 45 dB(A))

Medical device

Product measurements (LxWxH): 46 x 78 x 108 mm

Product weight: approx. 109 g

5 year guarantee

Sales unit: 3 / Shipping carton: 4 x 3

EAN no.: 4211125 60142 3

Item no.: 601.42