




BF 1000
Super Precision

beurer
wellbeing

BF 1000 Super Precision Diagnostic bathroom scale



Precise full-body analysis with 8 electrodes



Segmented display of body fat and muscle percentage
(arms, legs, trunk)

Cross measurement with 8 electrodes for precise monitoring of
different areas of the body (arms, legs, trunk)

Determination of body weight, body fat, body water, muscle
percentage, bone mass, BMR/AMR (kcal)

Determination of visceral fat (organic fat / abdominal fat)**

Basic turnover and active turnover in kcal

Raised dot matrix display on an aluminium stand with blue
backlight for optimum readability

Time and date display

Intuitive user guide in 5 languages (D/GB/F/ES/I)

With "beurer HealthManager Pro" / "beurer HealthManager" app
to download for free

With brushed stainless steel electrodes (handheld)

ITO*** coated weighing surface made from 10 mm thick safety glass

Change over to kg/lb/st

Capacity: 200 kg

Graduation:

- 0.1% for body fat, body water and muscle ratio

- 50 g for body weight and bone mass

5 activity levels

10 user memories with 30 measurement spaces each + guest mode

Overload indication

Button-On and Quick-Start technology

Automatic switch-off

Incl. 4 x 1.5V batteries AA

Product measurements (LxWxH): 38 x 39 x 100 cm

Product weight: approx. 6.35 kg

5 year guarantee

Sales unit: 1 / Shipping carton: -

EAN no.: 4211125 74860 9

Item no.: 748.60

* tested online application "beurer HealthManager" Cloud

** surrounds inner organs in the abdominal area. May lead to diseases like hypertension and
diabetes.

*** Indium Tin Oxide




748.60_BF1000_2022-01-17_02_DS_EN_BEU Subject to errors and changes




BF 1000
Super Precision

beurer
wellbeing

BF 1000 Super Precision Diagnostic bathroom scale



CONNECTION



Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via **Bluetooth®**



Syncs with
• iOS ≥ 12.0
• Android™ ≥ 8.0
• Bluetooth® 4.0



List of supported devices

748.60_BF1000_2022-01-17_02_DS_EN_BEU Subject to errors and changes

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.