


corefit I - EMS training belt



EMS technology



4 electrodes made of conductive carbon material



9 training programs



Stimulation of front and side abdominal muscles

4 electrodes (2x frontal abdominal muscles & 2x lateral abdominal muscles)

Easy to use

Flexible abdominal belt with hook and loop fastener

For waist sizes from 70 to max. 140 cm

Flexible ergonomic shape

LCD display with figurative symbols

9 training programs (20 to 30 min)

Intensity adjustable from 0 up to 40

Countdown timer with remaining time

Safety switch-off / pause function / automatic switch-off / key lock

Detachable control element

Flexible: 2 adhesive electrodes can be connected on the control for individual muscle stimulation (item no.: 647.05)

Battery changing indicator

The device must not be used by users of cardiac pacemakers, epileptics, pregnant women, people with cardiac arrhythmias

Output voltage: max. 70V p-p (500 Ohm load)

Output current: max. 140 mA p-p (500 Ohm load)

Output frequency: 5-85Hz

Pulse width: 200-400 μ s per phase

Pulseform: symmetrical, dual-phase square pulses

Incl. 3 x 1.5V batteries AAA

Product measurements control element (LxWxH): 8,8 x 9,1 x 2,2 cm

Product weight: approx. 330g (without batteries)

3 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 421112588105 4

Item no.: 881.05