



BM 53


beurer  
medical

## BM 53 Upperarm Blood pressure monitor

german  
engineering



Advanced arrhythmia detection:  
Detects atrial fibrillation for  
stroke prevention



Triple measurement is recommended by  
doctors for even more accurate  
measurements



Reliable measurement,  
even during pregnancy



Fully automatic blood pressure and pulse measurement on the upper arm

Patented resting indicator\* for exact measurement results

Cuff position check: For ensuring the cuff is correctly positioned on the upper arm

Arrhythmia detection: detects and warns of possible cardiac arrhythmias

Risk Indicator: Classification of measurements with coloured scale

Two user memories for 100 measured values each

Easy-to-read XL display with backlight

Calculation of the average value of all saved measured values and morning and evening blood pressure for the last 7 days

Automatic pressure preselection and deflation

Medical device

Date and time /automatic switch-off

Battery changing indicator

With premium cuff for upper arm circumferences from 22 - 42 cm

Incl. 4 x 1,5V AA batteries and storage pouch

Medical device according to new Medical Device Regulation (EU) 2017/745 (MDR)

Product measurements (LxWxH): 140 mm x 94 mm x 46 mm

Product weight: approx. 437 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65459 7

Item no.: 654.59

\* Patent EP 1 673 009

654.59\_BM53\_2023-10-09\_02\_DS\_EN\_BEU Subject to errors and changes