



BF 980 WIFI

beurer
wellbeing

BF 980 WIFI

Diagnostic bathroom scale WIFI / Bluetooth®



With app for downloading:
Analysis records and graphical
evaluations of body values



Incl. heart rate measurement



Integrated battery
for easy recharging



Quick & easy: *Bluetooth®* for initial set-up / WIFI for automatic data transfer

Automatic synchronisation: the data is automatically uploaded via WIFI each time you weigh yourself

Includes buttons for use without the app

ITO* coated weighing surface made of safety glass

Extra-large standing surface 35 x 30 cm

With extra large magic display in state-of-the-art dot matrix look

Digit height: 48 mm

Calculates body weight, body fat, muscle percentage and heart rate; other values via the app (metabolic age, body water, BMR/AMR calorie display, bone mineral mass)

5 activity levels

8 user memories with automatic user recognition

Change over kg/lb/st

Capacity: 180 kg

Graduation: 0.1 % for body fat and muscle percentage / 100 g for body weight / 1 bpm for heart rate

Overload indication

Automatic switch-on and switch-off function

Integrated battery for easy recharging (incl. USB-C cable)

Product measurements (LxWxH): 35 x 30 x 2.6 cm

Product weight: ca. 2.2 kg


5 year guarantee

Sales unit: 3 / Shipping carton: -

EAN No.: 4211125760090

Item No.: 760.09

*Indium tin oxide





760.09_BF980 WIFI_2022-01-27_01_DS_EN_BEU Subject to errors and changes



BF 980 WIFI

beurer
wellbeing

BF 980 WIFI
Diagnostic bathroom scale **WIFI / Bluetooth®**



Innovative networking between smartphone and scale - keep track of your measurements and body data at home and on the move. Easily synchronised with your smartphone via *Bluetooth®* / *WIFI*



• iOS ≥ 12.0
• Android™ ≥ 8.0
• Bluetooth® ≥ 4.2



List of supported devices
Liste der unterstützten Geräte
Liste des équipement pris en charge
Lista de los equipo compatibles
Elenco degli attrezzatura supportati

760.09_BF980 WIFI_2022-01-27_01_DS_EN_BEU Subject to errors and changes

The *Bluetooth®* word mark and logos are registered trademarks owned by *Bluetooth SIG, Inc.* and any use of such marks by *Beurer GmbH* is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of *Apple Inc.*, registered in the U.S. and other countries. *App Store* is a service mark of *Apple Inc.*, registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of *Google LLC*.

Android is a trademark of *Google LLC*.