




## FT 95 Non-contact thermometer


 german  
 engineering


Innovative networking between smartphone and thermometer – With free "beurer HealthManager" and "beurer BabyCare" app



With *Bluetooth*® low energy technology – Keep track of your measurements at home and on the move



LED temperature alarm (green/yellow/red)



Contactless infrared clinical thermometer

Measures without skin contact - hygienic and safe

Measurement of forehead, object and room temperature in seconds

Particularly suitable for taking the temperature of babies and small children

Measurement range:

Forehead temperature measurement: 34 °C bis 42.2 °C (93.2 °F bis 108 °F)

Object temperature measurement: 0 °C bis 80 °C (32 °F - 176 °F)

Measurement accuracy:

Forehead temperature measurement: ± 0.2 °C (± 0.4 °F)

Object temperature measurement: ± 4% oder ± 2 °C (± 4 °F)

60 memory spaces (forehead measuring)

°C/°F can be set

Display illumination (On/Off)

Sound control (On/Off)

Automatic switch off

Date and time

Incl. storage box

With battery changing indicator

Incl. 2 x 1.5V batteries AAA

Product measurements (LxWxH): 167.4 x 41.5 x 28 mm

Product weight: approx. 73.4 g

Medical device

5 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 4211125 79511 5

Item no.: 795.11



< 37,5°C  
everything ok




37,5-37,9°C  
increased temperature



≥ 38°C  
fever alert



## CONNECTION



\* tested online application "beurer HealthManager" Cloud

795.11-FT95-2018-08-22\_02\_D5\_EN\_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.