



## EM 37 Abdominal Muscle Belt

 german  
engineering



Electrical Abdominal Muscle  
Stimulation (EMS)



4 Electrodes with water contact  
(no contact gel, no wearing electrodes)



Stimulation of front and side abdominal muscles

Easy to use

Flexible abdominal belt with hook and loop fastener

Flexible stomach belt for waist sizes from 70 to max. 140 cm

Flexible ergonomic shape

2 adhesive electrodes incl. connecting cable for individual stimulation  
available as accessories Item no.: 647.05

LCD display with figurative symbols

5 training programs (22 to 31 min)

Intensity adjustable from 0 up to 40

Countdown timer with remaining time

Training memory function for continuous training control

Safety switch-off/ pause function/ automatic switch-off/ key lock

Removable control element

Belt hand washable

Battery changing indicator

The device must not be used by users of cardiac pacemakers, epileptics,  
pregnant woman, people with cardiac arrhythmias

Output voltage max. 70V p-p (500 Ohm load)

Output current: max. 140 mA p-p (500 Ohm load)

Output frequency: 30–80Hz

Pulse width: 200 µs per phase

Pulseform: symmetrical, dual-phase square pulses

Incl. 3 x 1.5V batteries AAA

Product measurements control element (LxWxH): 8.8 x 9.1 x 2.2 cm

3 year guarantee

Sales unit: 3 / Shipping carton: 6 x 3

EAN no.: 4211125 64703 2

Item no.: 647.03