



HA 50 Hearing amplifier


 german
 engineering

 barely
 visible


Comfortable fit behind the ear - barely visible

 low-noise
 amplification
 and rendering


Low-noise rendering and amplification of
 the volume of all sounds



Ideal for restricted hearing ability

Amplifies the volume of all sounds indoors and outdoors

Continuously adjustable volume

3 attachments to individually adjust to the ear canal

Incl. battery and replacement battery

Sound amplification: max. 40 dB

Maximum audio output: 128 dB

Frequency range: 100 - 6000 Hz

Medical product

Product measurements: 40 x 40 x 10 mm

Product weight: approx. 6 g (incl. battery and earplug)

Delivered in sales display

3 year guarantee

Sales unit: 8 / Shipping carton: 4 x 8

EAN no.: 4211125 64109 2

Item no.: 641.09



Included in delivery

Hearing amplifier
 3 earplugs in different sizes
 2 Typ 13/PR48 Batterien
 1 box for storage
 Instructions for use



Unit description

1. Microphone
 2. Volume control 1-4 (1 Min., 4 Max.)
 3. On/Off switch
 („I“ = on, „O“ = off)
 4. Headphone
 5. Earplug
 6. Battery compartment

