



BM 85

beurer
medical

BM 85 Bluetooth® - Upper arm blood pressure monitor



Resting indicator* for exact measurement result



With cuff for upper arm circumferences
from 22 - 42 cm



- Fully automatic blood pressure and pulse measurement on the upper arm
- Inflation Technology: Fast and convenient measurement, even during pumping up
- LED risk indicator: Coloured LED scale for classifying measurements
- Arrhythmia detection: Warning in case of possible heart rhythm disturbances
- 2 x 60 memory spaces
- Average value of all stored measurements and morning and evening blood pressure for the last 7 days
- Easy to use via illuminated sensor touch buttons
- White illuminated XL display
- Date and time
- Automatic switch-off
- Signal in case of application error
- Automatic pressure preselection and deflation
- Incl. USB power mains part, cuff holder and storage pouch
- With lithium-ion battery
- Medical device
- Mains operation:
- Input: 100-240 V~; 50-60 Hz; 60-120 mA
- Product measurements (LxWxH): 180 x 99 x 40 mm
- Product weight: approx. 456 g (without batteries, with cuff)
- 5 year guarantee
- Sales unit: 4 / Shipping carton: 2 x 4
- EAN no.: 4211125 65803 8
- Item no.: 658.03

*Patent EP 1 673 009



LED risk indicator

• Coloured LED scale for classifying measurements




Sensor touch

• Easy to use via illuminated sensor touch buttons




USB port

• USB interface for innovative networking between the blood pressure monitor and "beurer HealthManager" PC software and charging the lithium-ion battery



CONNECTION



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.



Syncs with
• iOS ≥ 10.0
• Android™ ≥ 4.4
• Bluetooth® ≥ 5.0



List of supported devices

658.03_BM85_2021-10-19_03_DS_EN_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

* tested online application "beurer HealthManager" Cloud
Can only be used with its full functionality with the „beurer HealthManager“ app.