



HA 20 Hearing amplifier


 german
 engineering


Amplifies the volume of all sounds indoors and outdoors



With Box for storage



Ideal for restricted hearing ability

Comfortable fit behind the ear

Continuously adjustable volume

3 attachments to individually adjust to the ear canal

Incl. battery and replacement battery

Sound amplification: max. 40 dB

Maximum audio output: 128 dB

Frequency range: 200 - 5000 Hz

Medical product

Product measurements: 40 x 40 x 12 mm

Product weight: approx. 9 g

Delivered in sales display

3 year guarantee

Sales unit: 8 / Shipping carton: 4 x 8

EAN No.: 4211125 64106 1

Item no.: 641.06



Included in delivery

- _____ Hearing amplifier
- _____ 3 earplugs in different sizes
- _____ 2 Typ 13/PR48 Batterien
- _____ 1 box for storage
- _____ Instructions for use



Unit description

- _____ 1. Microphone
- _____ 2. Volume control 1-4 (1 Min., 4 Max.)
- _____ 3. On/Off switch
 („I“ = on, „O“ = off)
- _____ 4. Headphone
- _____ 5. Earplug
- _____ 6. Battery compartment

