



BM 93  
CARDIO  
**beurer**  
medical

## BM 93 **CARDIO** Blood pressure monitor with ECG function

 german  
engineering

 **Bluetooth**<sup>®</sup>  
beurer  
CONNECT 



1-channel ECG device for recording  
heart rhythm



Advanced arrhythmia detection:  
Detects atrial fibrillation



With cuff for upper arm circumferences  
from 23 - 43 cm



Integrated ECG-electrodes

Optimum control of the blood pressure and heart rhythm via app

Fully automatic blood pressure and pulse measurement on the upper arm

Risk Indicator: Classification of measurements

Easy-to-read XL display

Two user memories for 60 measured values each

Alarm function

Medical device

Date and time /automatic switch-off

Battery changing indicator

Incl. 4 x 1,5V AAA batteries and storage pouch

Product measurements (LxWxH): Blood pressure monitor: 128 x 127 x 42 mm

Product weight: Blood pressure monitor: ca. 420 g (without batteries, with cuff)

5 year guarantee

Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125652296


Item no.: 652.29






BM 93  
CARDIO  
**beurer**  
medical




## CONNECTION




sync via

.....



Innovative networking between smartphone and blood pressure monitor - keep track of your measurements at home and on the move. Easily synchronised with your smartphone via Bluetooth®.

works with




HealthManager Pro

works with




HealthManager

works with




TÜV Rheinland CERTIFIED \*

Download on the



App Store


GET IT ON



Google Play

Syncs with

- iOS ≥ 12.0
- Android™ ≥ 8.0
- Bluetooth® ≥ 4.0



List of supported devices

652:29\_BM93\_2021-03-29\_02\_DS\_EN\_BEU Subject to errors and changes

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.

\* tested online application "beurer HealthManager" Cloud  
Can only be used with its full functionality with the "beurer HealthManager Pro" / "beurer HealthManager" app.